



Fluidnews

Pipettes | Automation | Pipette Service

Autumn / Winter 2017 - 2018



40% OFF PIPETMAN L RANGE

33% OFF ELECTRONIC PIPETMAN M

50% OFF RELOAD TIPS

15% OFF DISPENSAN INTRODUCTION OFFER

NEW CRISPR TECHNICAL NOTE

GRAND OPENING CELEBRATIONS



OUR NEW UK CENTRE OF EXCELLENCE OFFICIALLY OPENS

We are delighted to announce that we officially opened our new world-class, purpose built UK Centre of Excellence facility in Dunstable on the 21st September.

At an exclusive Grand Opening Champagne Reception, Andrew Selous, MP for South West Bedfordshire and Peter Tudor, Gilson's European Director, officially unveiled the new building.

It is our continued growth within the UK and our commitment to customers that has fuelled the demand and need to relocate to larger and more modern premises. Formerly housed at Charles Street in Luton for over thirty years, this new state of the art facility has been home to our amazing team of liquid handling experts, customer support and service technicians since we moved in April.



What's New?

We have invested in a new UKAS Accredited Pipette Service Centre which provides unrivalled support, calibration and maintenance of all makes of pipettes, not only Gilson's prestigious PIPETMAN brand.

Our new facility also features new application development laboratories and a larger, more streamlined distribution centre, enabling us to provide enhanced technical and customer support.

"This investment reinforces our commitment to science and provides the means necessary to support our continued growth within the UK and Europe.

Our mission is to partner with the scientific community and help researchers advance the pace of discovery. We do this by creating laboratory tools that improve the results of tomorrow and a solid infrastructure through our global Centres of Excellence that enables us to distribute our solutions and provide the best technical support and services possible.

Our new, modern facility compliments the recent evolution of our corporate branding, which is aligned to the company's path and vision to enable verifiable science and make lab life easier for researchers."

Peter Tudor, Gilson's European Director

Champagne Reception

The event was celebrated with key customers, partners and members of the business community enjoying a champagne reception, lunch and speeches

A big thank you to all our customers who have helped us develop Gilson into the industry leading organisation that it is today!





pipetman®

L

INCREASE YOUR PRODUCTIVITY

All Day Pipetting Comfort

- Ideal weight, handle shape and tip ejection for maximum comfort
- Pipette for longer without fatiguing your hand
- Volume-locking mechanism prevents accidental volume change
- Ideal for even the most intensive pipetting routines
- Now **fully autoclavable without disassembling**



40% OFF
PIPETMAN L RANGE
& STARTER KITS
 OFFER CODE: GIL133D



40% OFF PIPETMAN L SINGLE CHANNELS			
MODEL / VOLUMES	PART NUMBERS		40% OFF
	S/S Ejector	Plastic Ejector	
P2L (0.2-2µL)	FA10001M	FA10001P	
P10L (1-10µL)	FA10002M	FA10002P	
P20L (2-20µL)	FA10003M	FA10003P	
P100L (10-100µL)	FA10004M	FA10004P	£221.00 £132.60
P200L (20-200µL)	FA10005M	FA10005P	
P1000L (100-1000µL)	FA10006M	FA10006P	
P5000L (500-5000µL)*		FA10007	
P10mL L (1-10mL)*		FA10008	
40% OFF PIPETMAN L MULTICHANNELS			
L8 x 10 (0.5-10µL)	FA10013		
L8 x 20 (2-20µL)	FA10009		£633.00 £379.80
L8 x 200 (20-200µL)	FA10011		
L8 x 300 (20-300µL)	FA10015		
L12 x 10 (0.5-10µL)	FA10014		
L12 x 20 (2-20µL)	FA10010		£739.00
L12 x 200 (20-200µL)	FA10012		£443.40
L12 x 300 (20-300µL)	FA10016		
40% OFF PIPETMAN L STARTER KITS			
Microvolume Kit P2L, P10L & P100L	F167550		£631.00 £378.60
PIPETMAN L Kit P20L, P200L & P1000L	F167350		£631.00 £378.60
PIPETMAN L 4-Pipette Kit P2L, P20L, P200L & P1000L	F167370		£779.00 £467.40

ORDER NOW: QUOTE OFFER CODE GIL133D

Call: 0845 5196 872 or email: orders@gilsonuk.com
 Order Online: www.gilsonuk.com/promotions

pipetman®

M

GUARANTEED PERFORMANCE & ACCURACY

Pipette With More Consistency & Less Effort

- Electronic pipette offering easy-to-use features
- Guarantees accuracy and precision in both standard and repetitive pipetting modes
- Lightweight and comfortable
- Virtually zero pipetting forces
- Streamlined user interface with five pipetting modes, including protocol customisation
- 33% off both single and multichannels



33% OFF
ELECTRONIC
PIPETMAN M RANGE
 OFFER CODE: GIL142D

33% OFF PIPETMAN M SINGLE CHANNELS			
MODEL / VOLUMES	PART NUMBERS	LIST PRICE	33% OFF
P20M (2-20µL)	F81023		
P100M (5-100µL)	F81013		
P200M (20-200µL)	F81024	£383.00	£256.61
P300M (20-300µL)	F81014		
P1200M (100-1200µL)	F81015		
P5000M (500-5000µL)	F81016		
P10mLM (1-10mL)	F81017		
33% OFF PIPETMAN M MULTICHANNELS			
P8 x 10M (0.5-10µL)	F81025		
P8 x 20M (1-20µL)	F81027		
P8 x 100M (10-100µL)	F81018	£758.00	£507.86
P8 x 200M (20-200µL)	F81029		
P8 x 300M (10-300µL)	F81031		
P8 x 1200M (50-1200µL)	F81020		
P12 x 10M (0.5-10µL)	F81026		
P12 x 20M (1-20µL)	F81028		
P12 x 100M (10-100µL)	F81019	£875.00	£586.25
P12 x 200M (20-200µL)	F81030		
P12 x 300M (10-300µL)	F81032		
P12 x 1200M (50-1200µL)	F81021		

ORDER NOW: QUOTE OFFER CODE GIL142D

Call: 0845 5196 872 or email: orders@gilsonuk.com
 Order Online: www.gilsonuk.com/promotions

PIPETMAN L and PIPETMAN M offers may not be used in conjunction with any other discount. Valid until 31/01/18 to UK mainland customers only.





pipetman[®]

DIAMOND TIPS

HIGHEST QUALITY, LOWEST PRICE



Save Space, Time and Money

- RELOAD PACKs are the easiest and quickest way to refill an empty TIPACK Box
- 3x smaller than a regular tip box and take up 60% less bench space for easy storage
- Eco-friendly, uses 70% less plastic than other racked boxes
- Only PIPETMAN tips guarantee the performance of your PIPETMAN pipettes for reliable results
- RELOAD Packs are one of the most cost effective ways of purchasing PIPETMAN tips

50% OFF
RELOAD TIPS
 OFFER CODE: GIL152D

50% OFF RELOAD PACKS			
MODEL / VOLUMES	PART NUMBERS	LIST PRICE	50% OFF
DL10 0.1-20µL (960 tips)	F167012	£27.00	£13.50
D200 2-200µL (960 tips)	F167013	£24.00	£12.00
D1000 100-1000µL (576 tips)	F167014	£16.00	£8.00

ORDER NOW: QUOTE OFFER CODE GIL152D

Call: 0845 5196 872 or email: orders@gilsonuk.com

Order Online: www.gilsonuk.com/promotions

dispensman[™]

NEW

DISPENSE LIQUIDS SIMPLY & ACCURATELY

Precise

Dispense just the liquid needed every time with a quick, easy, and precise volume setting.

No Loss of Liquid

Prime your operation by purging the air in RECYCLE position keeping the liquid inside the bottle and avoiding any loss of liquid.

Safe

Keep your bench and other devices safe with the anti-drip system and the return of liquid into the bottle after use.

Flexible

- Numerous bottle adaptors
- Adaptable to every type of liquid with an easy calibration by the user
- Fully autoclavable at 121°C
- Convenient 3-position nozzle



15% OFF INTRODUCTORY OFFER DISPENSAN				
MODEL	PART NUMBER	VOLUME RANGE	SCALING (ML)	15% OFF
DISPENSAN 2.5 mL	F110101	0.25-2.5mL	0.05	£220.00
DISPENSAN 5.0 mL	F110102	0.5-5mL	0.10	£187.00
DISPENSAN 10 mL	F110103	1-10mL	0.20	£245.00
DISPENSAN 25 mL	F110104	2.5-25mL	0.50	£245.00
DISPENSAN 50 mL	F110105	5-50mL	1.00	£208.25

15% OFF
INTRODUCTORY OFFER
 OFFER CODE: GIL153D

ORDER NOW QUOTE OFFER CODE GIL153D

Call: 0845 5196 872 or email: orders@gilsonuk.com

Visit: www.gilsonuk.com/products/pipettes/large-volume.html





platemaster®

FAST MANUAL MICROPLATE PIPETTING



PLATEMASTER is an amazingly easy-to-use, compact and accurate solution for high throughput pipetting of 96 and 384 well microplates.

Save Time And Increase Your Productivity

With minimal setup you can fill 96 or 384 well plates in less than one minute by reducing the number of manipulations for filling a microplate.

Simplifies Sample Prep

PLATEMASTER requires no programming, no power and minimal training, so you can start pipetting right away.

Convenient For Any Microplate Application

Achieve accuracy and precision on the bench or under the hood.

ORDER NOW OR ARRANGE A DEMO

Call: 0845 5196 872 or email: orders@gilsonuk.com

Trade-In & Upgrade

REPLACE YOUR OLD LIQUID HANDLER OR HPLC SYSTEM



Trade-in any make and model of automated liquid handler or HPLC system and get 25% off a new Gilson instrument.

25% Off Automated Liquid Handlers

- Gilson offer small footprint to larger scale solutions
- Flexible and versatile platforms with dependable functionality and performance
- Task orientated TRILUTION® software and customisable bed layouts

25% Off PLC Purification Systems

- HPLC and FLASH all in one instrument
- Compact design and cost effective
- Three models available -flow rates up to 500mL/min
- Customisable with advanced detection options

FIND OUT MORE OR REQUEST QUOTE

Call: 0845 5196 872 or email: info@gilsonuk.com and we will help you trade-in your old system.

Visit gilsonuk.com and search "platemaster"

£7,500

with free tips to the value of £500

OFFER CODE: CP0362

Valid until 31/01/18 to UK customers only.

25% OFF

when you trade-in
Email our team on
info@gilsonuk.com

Valid until 31/01/18 to UK customers only.






pipetmax®

MEET MAX YOUR NEW LABORATORY ASSISTANT

Do you suffer from Repetitive Strain Injury (RSI)?

Do you perform the same repetitive assay, day after day?

LET PIPETMAX TAKE THE STRAIN!



PIPETMAN Inside

Search "pipetmax" on gilsonuk.com



Liquid Handling from the Experts

With a PIPETMAN® inside, you can be assured of world class pipetting, giving you the kind of accuracy and precision you expect. With volumes from 0.5µL to 1mL, PIPETMAX can perform liquid transfers from tubes, vials, plates and reservoirs, and can handle both 96 and 384 well plates.



PIPETMAN® is an easy-to-use, automated liquid handling platform that performs the mundane task of sample preparation for you! Enjoy more freedom away from the bench to focus on what's important; value-added tasks and your next discovery.

With a multitude of applications at your disposal, PIPETMAX increases the reproducibility and pace of your experiments.

Suitable for ALL Life Science Applications

From qPCR and CRISPR to ELISA and cell-based assays, PIPETMAX is suitable for a wide variety of applications. Software wizards for qPCR and Normalization make set up easy and use a language spoken by scientists. Magnetic bead racks and heating and shaking options, expand the repertoire of PIPETMAX over and above other liquid handlers.

GET YOUR OWN MAX MASCOT

Email our team on info@gilsonuk.com or call 0845 5196 872

FIND OUT MORE OR REQUEST A QUOTE

Call: 0845 5196 872 or email: info@gilsonuk.com to arrange a demo or get a quote

NEW CRISPR TECHNICAL NOTE

CRISPR/Cas9 technology is a cutting-edge genetic-engineering tool, which allows the rapid and precise editing of genetic sequences in a wide range of organisms.

REQUEST A COPY OF THE GILSON CRISPR TECHNICAL NOTE

Get your copy of the Gilson TECHNICAL Note: *"CRISPR/Cas9 Genetic Engineering - How Gilson Manual and Automated Pipetting Tools Facilitate the Process"* to learn more and see how Gilson are providing innovative products to enable pioneering new science.

Email: info@gilsonuk.com



GILSON PIPETTE SERVICE CENTRE

Our new UKAS accredited, purpose-built pipette service centre offers unrivalled service, maintenance and calibration for all makes and models of pipettes.

- Fully trained technical professionals who are pipetting experts
- We are UKAS accredited and offer UKAS calibration
- Onsite InLab™ Pipette Service Clinics by our qualified technicians
- Rapid turnaround for Back-to-Base returned pipettes
- Drop Box pipette service scheme at selected locations
- Trained to service all makes and models of pipettes

1-DAY PIPETTING TRAINING COURSE

Join the Gilson family of experts and become a proficient and skilled certified pipetting user. With a series of 1-day educational training courses, learn essential skills and knowledge that will allow you to pipette accurately and effectively. Topics include:

- Good pipetting practice and verification procedures
- Manual and automated pipetting techniques
- Basic pipette verification & maintenance training
- How to pipette difficult liquids
- Reverse and touch off pipetting



ARRANGE A SERVICE

- Call: 0845 643 5216
- E: servicesales@gilsonuk.com
- Visit: www.gilsonuk.com

LIKE TO ATTEND?

Register your interest by:

- Calling: 0845 5196 872
- Emailing: info@gilsonuk.com

Gilson UK
3B Humphrys Road
Woodside Estate
Dunstable
LU5 4TP

www.gilsonuk.com
info@gilsonuk.com